
本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human RSPO3 Protein (aa 1-146, His Tag)

产品货号: TP06224

第三版

别名:CRISTIN1;PWTSR;THSD2

描述:R-spondin 3 (RSPO3) is a member of the R-Spondin (RSPO) family in vertebrates that activate Wnt/beta-catenin signaling, plays a key role in these processes. The RSPO family of secreted Wnt modulators is involved in development and disease and holds therapeutic promise as stem cell growth factors. The four members have high structural homology. RSPO2 and RSPO3 are more potent than RSPO1, whereas RSPO4 is relatively inactive. All RSPO members require Wnt ligands and LRP6 for activity and amplify signaling of Wnt3A, Wnt1, and Wnt7A, suggesting that RSPO proteins are general regulators of canonical Wnt signaling. RSPO3/PCP signaling during gastrulation requires Wnt5a and is transduced via Fz7, Dvl, and JNK. RSPO3 functions by inducing Sdc4-dependent, clathrin-mediated endocytosis. RSPO3 is a novel, evolutionarily conserved angiogenic factor in embryogenesis. RSPO3 has a key role in the interaction between chorion and allantois in labyrinthine development.

配方:PBS

分子量:15.3 kDa

序列:Met 1-Val 146

纯度:> 95% by HPLC

浓度:

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.