

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human MAGOH Protein (His tag)

产品货号: TP06289

第三版

别名:MAGOH1;MAGOHA

描述:Drosophila that have mutations in their mago nashi (grandchildless) gene produce progeny with defects in germlasm assembly and germline development. This gene encodes the mammalian mago nashi homolog. In mammals, mRNA expression is not limited to the germ plasm, but is expressed ubiquitously in adult tissues and can be induced by serum stimulation of quiescent fibroblasts.

配方:PBS

分子量:18 kDa

序列:Met 1-Ile146

纯度:> 95% by HPLC

浓度:

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.