

---

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

---

Rat IL-22R / IL22RA1 Protein (ECD, Fc Tag)

产品货号: TP06353

第三版

**别名:**IL-22R, IL22RA1

**描述:**IL-22R belongs to the type II cytokine receptor family. It contains 2 fibronectin type-III domains and is expressed in colon, liver, lung, pancreas and kidney. IL-22R also can be expressed in keratinocytes of normal skin as well as in psoriatic skin lesion. Overexpression of IL-22R can be detected in synovial fluid cells from rheumatoid arthritis and spondyloarthritis patients. IL-22R is a component of the receptor for IL20, IL22 and IL24. The component of IL-22R formed by IL22RA1 and IL10RB enables IL22 signaling via JAK/STAT pathways. IL22 also induces activation of MAPK1/MAPK3 and Akt kinases pathways. Component of one of the receptor for IL20 and IL24 formed by IL22RA1 and IL20RB also signaling through STATs activation. IL-22R mediates IL24 antiangiogenic activity as well as IL24 inhibitory effect on endothelial cell tube formation and differentiation.

**配方:**PBS**分子量:**50.8 kDa**序列:**Met1-Ala228**纯度:**> 95% by HPLC**浓度:****内毒素:**<1.0 EU per 1 ug of protein (determined by LAL method)**存储:**+4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.