

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human FLT3L / Flt3 ligand / FLT3LG Protein (His tag)

产品货号: TP06360

第三版

别名:FL;FLT3L;FLT3LG

描述:FLT3L, also known as flt3 ligand, is a small molecule that acts as a growth factor that increases the number of immune cells by activating the hematopoietic progenitors. In vivo, FLT3L also induces the mobilization of the hematopoietic progenitors and stem cells. This may help the system to kill cancer cells. Dendritic cells (DCs) provide the key link between innate and adaptive immunity by recognizing pathogens and priming pathogen-specific immune responses. FLT3L controls the development of DCs and is particularly important for plasmacytoid DCs and CD8⁺-positive classical DCs and their CD103⁺-positive tissue counterparts.

配方:PBS

分子量:20.2 kDa

序列:Thr 27-Pro185

纯度:> 95% by HPLC

浓度:

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.