

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Rat HepaCAM2 Protein (Fc tag)

产品货号: TP06376

第三版

别名:HepaCAM2

描述: HepaCAM2 belongs to the immunoglobulin family. HepaCAM2, together with HEPACAM1, are cell adhesion molecules. HEPACAM1 modulates cell adhesion and migration. It also inhibits cancer cell growth. HepaCAM2 contains 2 Ig-like C2-type (immunoglobulin-like) domains. The immunoglobulin domain is a type of protein domain that consists of a 2-layer sandwich of between 7 and 9 antiparallel β -strands arranged in two β -sheets with a Greek key topology. Members of the immunoglobulin superfamily are found in hundreds of proteins of different functions.

配方: PBS

分子量: 62.3 kDa

序列: Met1-Leu350

纯度: > 95% by HPLC

浓度:

内毒素: <1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.