
本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human CNPY4 Protein (His tag)

产品货号: TP06386

第三版

别名: PRAT4B

描述: CNPY4 belongs to the canopy family. CNPY4 interacts with toll-like receptor 4 (TLR4) and plays a role in the regulation of the cell surface expression of TLR4. Toll-like receptors (TLRs) recognize microbial products and induce immune responses. Lipopolysaccharide is recognized by the receptor complex consisting of TLR4 and MD-2. As CNPY4, PRAT4B also regulates cell surface expression of TLR4. PRAT4B has a signal peptide followed by a mature peptide. It is associated with the hypoglycosylated, immature form of TLR4 but not with MD-2 or TLR2.

配方: PBS

分子量: 27.4 kDa

序列: Met 1-Leu248

纯度: > 95% by HPLC

浓度:

内毒素: <1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.