

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human BCL6 / B-cell CLL lymphoma 6 Protein (aa 1-150, His & Trx tag)

产品货号: TP06399

第三版

别名: BCL5; BCL6A; LAZ3; ZBTB27; ZNF51

描述: The protein encoded by this gene is an evolutionarily conserved 95-kDa protein containing six C-terminal zinc-finger motifs and an N-terminal POZ domain. It has been reported that BCL-6 is present in DNA-binding complexes in nuclear extracts from various B-cell lines. There are many relationships between non-Hodgkin's lymphoma, diffuse large cell lymphoma and BCL6 's translocations. BCL6 can repress transcription from promoters linked to its DNA target sequence and this activity is dependent upon specific DNA-binding and the presence of an intact N-terminal half of the protein.

配方: PBS

分子量: 34.2 kDa

序列: Met 1-Met 150

纯度: > 95% by HPLC

浓度:

内毒素: <1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.