

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human AIM2 Protein (GST tag)

产品货号: TP06401

第三版

别名:PYHIN4

描述:AIM2, Absent In Melanoma 2 is a member of the interferon-inducible HIN-200 protein family that contains an amino-terminal pyrin domain and a carboxy-terminal oligonucleotide/oligosaccharide-binding domain, senses cytoplasmic DNA by means of its oligonucleotide/oligosaccharide-binding domain and interacts with ASC (apoptosis-associated speck-like protein containing a CARD) through its pyrin domain to activate caspase-1. In response to foreign cytoplasmic DNA, AIM2 forms an inflammasome, resulting in caspase activation in inflammatory cells. It had been pointed to a role of AIM2 function in both inflammation and cancer. AIM-2 antigen is expressed in a wide variety of tumor types, including neuroectodermal tumors, as well as breast, ovarian and colon carcinomas. AIM-2 could be used as a tumor antigen target for monitoring vaccine trials or to develop antigen specific active immunotherapy for glioma patients.

配方:PBS

分子量:65.2 kDa

序列:Met 1-Thr 343

纯度:> 95% by HPLC

浓度:

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.