

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human SHISA3 Protein (Fc Tag)

产品货号: TP06403

第三版

别名:hShisa3

描述:SHISA3 encodes a single-transmembrane protein which is one of nine members of a family of transmembrane adaptors that modulate both WNT and FGF signaling by blocking the maturation and transport of their receptors to the cell surface. Members of this family contain an N-terminal cysteine-rich domain with a distinct cysteine pattern, a single transmembrane domain, and a C-terminal proline-rich, low complexity region. The encoded protein acts as a tumor suppressor by accelerating beta-catenin degradation.

配方:PBS

分子量:34.9 kDa

序列:Met 1-Pro 98

纯度:> 95% by HPLC

浓度:

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.