

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Mouse Meteorin / METRN Protein (Fc tag)

产品货号: TP06454

第三版

别名:1810034B16Rik;Hyrac

描述:Meteorin is a novel secreted protein that is expressed in undifferentiated neural progenitors and in the astrocyte lineage, including radial glia. It plays important roles in the differentiation of glial cells and also in axonal network formation during neurogenesis. Meteorin selectively promoted astrocyte formation from mouse cerebrocortical neurospheres in differentiation culture, whereas it induced cerebellar astrocytes to become radial glia. Meteorin also induced axonal extension in small and intermediate neurons of sensory ganglia by activating nearby satellite glia.

配方:PBS

分子量:58.4 kDa

序列:Gly22-Asp291

纯度:> 95% by HPLC

浓度:

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.