



本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human SLITRK4 / Slitrk4 Protein (His tag)

产品货号: TP06468

第三版

别名:DKFZp547M2010

描述:SLIT and NTRK-like family proteins 4 is a member of the SLIT and NTRK-like family. This protein is characterized as an integral membrane protein that has two leucine-rich repeat (LRR) domains and a carboxy-terminal domain that is partially similar to SLIT protein and TRK Neurotrophin. There are six members in this family: SLITRK1, SLITRK2, SLITRK3, SLITRK4, SLITRK5, and SLITRK6. Their expression was detected mainly in the brain, but the expression profile of each SLITRK was unique. Each of these members is located in nervous system except for SLITRK6. SLITRK expression was also investigated in various types of brain tumor tissue.

配方:PBS

分子量:68.6 kDa

序列:Met 1-Pro 616

纯度:> 95% by HPLC

浓度:

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.