

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Mouse CREG / CREG1 Protein (Fc tag)

产品货号: TP06473

第三版

别名:AA755314;Creg

描述:CREG1 belongs to the CREG family. It is a adenovirus E1A protein which both activates and represses gene expression to promote cellular proliferation and inhibit differentiation. Thus it may contribute to the transcriptional control of cell growth and differentiation. The transcriptional control activity of cell growth requires interaction with IGF2R. CREG1 also antagonizes transcriptional activation and cellular transformation by E1A. It shares limited sequence similarity with E1A and binds both the general transcription factor TBP and the tumor suppressor pRb in vitro. CREG1 gene may contribute to the transcriptional control of cell growth and differentiation.

配方:PBS

分子量:48.5 kDa

序列:Met1-Gln220

纯度:> 95% by HPLC

浓度:

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.