



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Rat IL9R / Interleukin 9 receptor Protein (His tag)

产品货号: TP06482

第三版

别名:IL9R

描述:Interleukin 9 (IL9) is involved in mast cell maturation and the enhancement of IgE production by B cells. Furthermore, linkage data in human and mice have suggested that IL9 may contribute to asthma. IL9R is known to be autosomal in mice, and is X-linked only in primates. The more recent X linkage and more telomeric position of the IL9R gene may explain its autosomal, 'uninactivated' transcriptional status.

配方:PBS

分子量:28.2 kDa

序列:Met 1-Ala 270

纯度:> 95% by HPLC

浓度:

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.