

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Mouse MFI2 / CD228 / melanotransferrin Protein (Fc tag)

产品货号: TP06487

第三版

别名:CD228;MTf

描述:The long non-coding RNA MFI2 antisense RNA is overexpressed in human cancer tissues and its increased expression is associated with occurrence and metastasis of cancer. lncRNA MFI2 could promote proliferation and migration of osteosarcoma cells by regulating FOXP4 expression, which suggested critical roles of lncRNA MFI2 and FOXP4 in occurrence and development of human osteosarcoma.

配方:PBS

分子量:103 kDa

序列:Met 1-Gln 708

纯度:> 95% by HPLC

浓度:

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C 或者-70 ° C. 避免反复冻融.