

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human PTPN12 Protein (aa 1-355, His & GST tag)

产品货号: TP06489

第三版

别名:PTP-PEST;PTPG1

描述:PTPN12 is a member of the protein tyrosine phosphatase (PTP) family. PTPs are known to be signaling molecules that regulate a variety of cellular processes including cell growth, differentiation, mitotic cycle, and oncogenic transformation. PTPN12 contains a C-terminal PEST motif, which serves as a protein – protein interaction domain, and may be related to protein intracellular half-life. PTPN12 was found to bind and dephosphorylate the product of oncogene c-ABL, thus may play a role in oncogenesis. PTPN12 was shown to interact with, and dephosphorylate, various of cytoskeleton and cell adhesion molecules, such as p130 (Cas), CAKbeta/PTK2B, PSTPIP1, and paxillin, which suggested its regulatory roles in controlling cell shape and mobility.

配方:PBS

分子量:69.4 kDa

序列:Met 1-Gln355

纯度:> 95% by HPLC

浓度:

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.