



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human CDC42BPB Protein (His & GST tag)

产品货号: TP06502

第三版

别名:MRCKB

**描述:**CDC42BPB is a member of the serine / threonine protein kinase family that contains a Cdc42 / Rsc-binding p21 binding domain similar to that of PAK kinase. The kinase domain of this protein is related to the myotonic dystrophy kinase related ROK and this kinase may have functions in downstream regulating of Cdc42 in cytoskeletal recognition. It has been reported that the CDC42BPB protein take part in regulating numerous cellular functions by binding to members of a serine / threonine protein kinase subfamily. These functions include the remodeling of the cell cytoskeleton that is a feature of cell growth and differentiation.

**配方:**PBS

**分子量:**82.4 kDa

**序列:**Met 1-His 427

**纯度:**> 95% by HPLC

**浓度:**

**内毒素:**<1.0 EU per 1 ug of protein (determined by LAL method)

**存储:** +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.