
本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human PTMA Protein (GST tag)

产品货号: TP06524

第三版

别名:TMSA

描述:PTMA (prothymosin, alpha, N-GST chimera) is a small, 12.4 kDa protein. It is a 109-111 amino acid long polypeptide as the precursor of thymosin a1. Thymosins are named because they were originally isolated from the thymus. But now in many other tissues, thymosins also can be detected. Thymosins have diverse biological activities, and two in particular, thymosins a1 and _4, have potentially important uses in medicine, some of which have already progressed from the laboratory to the clinic. In general, PTMA is associated with cellular proliferation and carcinogenesis (Eschenfeldt et al., 1986), cellular and viral transcription (Cotter et al., 2000), protection against apoptosis and chromatin remodelling (Karetsou et al., 1998). PTMA may have a dual role both intracellularly and extracellularly. In relation to diseases, thymosins have been categorized as biological response modifiers. Thymosin a1 is derived from PTMA. For animals that lack thymus glands, thymosin a1 is responsible for the activity of that preparation in restoring immune function.

配方:PBS

分子量:39 kDa

序列:Ser 2-Asp 111

纯度:> 95% by HPLC

浓度:

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.