
本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Mouse CSRP1 / CSRP / CRP1 Protein (His tag)

产品货号: TP06528

第三版

别名:AA408841;AA959891;AW545626;CRP1;Csrp

描述:Cysteine and glycine-rich protein 1, also known as Cysteine-rich protein 1, CSRP1 and CSRP, is a member of the CSRP family which may be involved in regulatory processes important for development and cellular differentiation. CSRP1 contains two LIM zinc-binding domains. The LIM / double zinc-finger motif found in CSRP1 is found in a group of proteins with critical functions in gene regulation, cell growth, and somatic differentiation. Zebrafish CSRP1 is expressed in the mesendoderm and its derivatives. CSRP1 interacts with Dishevelled 2 (Dvl2) and Diversin (Div), which control cell morphology and other dynamic cell behaviors via the noncanonical Wnt and JNK pathways. When CSRP1 message is knocked down, abnormal convergent extension cell movement is induced, resulting in severe deformities in midline structures. In addition, cardiac bifida is induced as a consequence of defects in cardiac mesoderm cell migration. CSRP1 acts as a key molecule of the noncanonical Wnt pathway, which orchestrates cell behaviors during dynamic morphogenetic movements of tissues and organs.

配方:PBS

分子量:22 kDa

序列:Met 1-Glu 193

纯度:> 95% by HPLC

浓度:

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.