

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Mouse TALLA-1 / TSPAN7 Protein (His tag)

产品货号: TP06532

第三版

别名:1200014P11Rik;A15;AI323365;Cd231;Mxs1;PE31;R74651;TALLA;Tm4sf2

描述:TALLA-1, also known as TSPAN7, is a member of the transmembrane 4 superfamily. Most members of this family are cell-surface proteins that are characterized by the presence of four hydrophobic domains. TALLA-1 gene is associated with X-linked mental retardation and neuropsychiatric diseases such as Huntington's chorea, fragile X syndrome and myotonic dystrophy. TALLA-1 is a cell surface glycoprotein and may have a role in the control of neurite outgrowth. It is known to complex with integrins.

配方:PBS

分子量:13 kDa

序列:Arg113-Met213

纯度:> 95% by HPLC

浓度:

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.