
本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human Resistin / ADSF / RETN Protein (Fc Tag)

产品货号: TP06533

第三版

别名:ADSF;FIZZ3;RETN1;RSTN;XCP1

描述:cGMP-dependent 3',5'-cyclic phosphodiesterase, also known as cyclic GMP-stimulated phosphodiesterase and PDE2A, is a peripheral membrane protein which belongs to the cyclic nucleotide phosphodiesterase family and PDE2 subfamily.

Phosphodiesterases (PDEs) comprise a family of enzymes that regulate the levels of cyclic nucleotides, key second messengers that mediate a diverse array of functions. Phosphodiesterases (PDEs) modulate signaling by cyclic nucleotides in diverse processes such as cardiac contractility, platelet aggregation, lipolysis, glycogenolysis, and smooth muscle contraction. PDE2A is an evolutionarily conserved cGMP-stimulated cAMP and cGMP PDE. PDE2A contains two GAF domains. PDE2A is expressed in brain and to a lesser extent in heart, placenta, lung, skeletal muscle, kidney and pancreas. PDE2A is a cyclic nucleotide phosphodiesterase with a dual-specificity for the second messengers cAMP and cGMP, which are key regulators of many important physiological processes. PDE2A is involved in the regulation of blood pressure and fluid homeostasis by the atrial natriuretic peptide (ANP), making PDE2-type enzymes important targets for drug discovery.

配方:PBS

分子量:38.2 kDa

序列:Ser17-Pro108

纯度:> 95% by HPLC

浓度:

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.