

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human Peroxiredoxin 6 / PRDX6 Protein (His tag)

产品货号: TP06558

第三版

别名:1-Cys;aiPLA2;AOP2;HEL-S-128m;NSGPx;p29;PRX

描述:acid phosphatase-like protein 2, also known as ACPL2, is a secreted protein which belongs to the histidine acid phosphatase family. A large-scale effort, termed the Secreted Protein Discovery Initiative (SPDI), was undertaken to identify novel secreted and transmembrane proteins. In the first of several approaches, a biological signal sequence trap in yeast cells was utilized to identify cDNA clones encoding putative secreted proteins. A second strategy utilized various algorithms that recognize features such as the hydrophobic properties of signal sequences to identify putative proteins encoded by expressed sequence tags (ESTs) from human cDNA libraries. A third approach surveyed ESTs for protein sequence similarity to a set of known receptors and their ligands with the BLAST algorithm. Finally, both signal-sequence prediction algorithms and BLAST were used to identify single exons of potential genes from within human genomic sequence.

配方:PBS

分子量:26.5 kDa

序列:Met 1-Pro 224

纯度:> 95% by HPLC

浓度:

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.