

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Mouse Carbonic Anhydrase X / Car10 Protein (His tag)

产品货号: TP06567

第三版

别名:2700029L05Rik;BB085816;Ca10;Car10;RP23-448I9.1

描述:Carbonic anhydrase X, also called carbonic anhydrase - related protein X (CARPX) and CA10, belongs to the CA family of zinc metalloenzymes which catalyze the reversible hydration of carbon dioxide in various biological processes such as respiration, renal tubular acidification and bone resorption. The secreted protein CARPX without CA activity (hydration of CO₂) is identified as an acatalytic member of the alpha-carbonic anhydrase subgroup. CARP X expression is detected in the adult total brain and almost all parts of the central nervous system, but not in the fetal brain. Accordingly, CARP X is suggested to play a role in the development of central nervous system, especially the brain. The same CARP X protein are encoded by multiple transcript variants of this gene.

配方:PBS

分子量:33.1 kDa

序列:Met 1-Asn 300

纯度:> 95% by HPLC

浓度:

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.