

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human NRN1 Protein (His Tag)

产品货号: TP06593

第三版

别名:dJ380B8.2;MGC44811;NRN

描述:Neuritin 1 (NRN1) is a member of neuritin family. Neuritin is a glycosylphosphatidylinositol- anchored protein induced by neural activity. It is expressed in postmitotic-differentiating neurons of the developing nervous system and a population of small-diameter neurons in the dorsal root ganglia and was anterogradely and retrogradely transported. Neuritin message is induced by neuronal activity and by the activity-regulated neurotrophins BDNF, nerve growth factor (NGF) and NT-3. Purified recombinant neuritin promotes neurite outgrowth and arborization in primary embryonic hippocampal and cortical cultures. Thus, neuritin is considered as a downstream effector of activity-induced neurite outgrowth. In clinical, neuritin levels in diabetes were reduced in both dorsal root ganglia and sciatic nerve of rats, and these deficits were reversed in vivo by treatment with NGF. This manipulation of neuritin levels in diabetes may provide a potential target for the therapeutic intervention in the management of neuropathy.

配方:PBS

分子量:

序列:Met 1-Asn 115

纯度:> 95% by HPLC

浓度:

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.