
本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human PS6K / RPS6KB1 Protein (GST tag)

产品货号: TP06605

第三版

别名:p70 S6KA;p70(S6K)-alpha;p70-alpha;p70-S6K;PS6K;RPS6KB1;S6K;S6K-beta-1;S6K1;STK14A

描述:PS6K, also known as RPS6KB1, is a serine/threonine-protein kinase. It belongs to the RSK (ribosomal s6 kinase) family. Members of this family function in signal transduction. PS6K is an isoform of p70 ribosomal S6 kinase (S6K). S6K can be activated by mitogenic stimuli such as growth factors, insulin and cytokines. It phosphorylates the ribosomal protein S6. PS6K also phosphorylates other proteins such as eIF4B, eEF2K and SKAR. It is a crucial effector of mTOR(rapamycin) signaling. PS6K is dissociated from the EIF3 complex and activated upon mitogenic stimulation, phosphorylation by the mammalian target of mTOR complex 1 (mTORC1). Its active form then phosphorylates and activates several substrates in the preinitiation complex, including the EIF2B complex and the cap-binding complex component EIF4B. PS6K also functions in cell proliferation, cell growth and cell cycle progression.

配方:PBS

分子量:85.4 kDa

序列:Met 1-Leu525

纯度:> 95% by HPLC

浓度:

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.