

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human Ube2L6 Protein (His Tag)

产品货号: TP06607

第三版

别名:RIG-B;UBCH8

**描述:**UBCH8, also known as UBE2L6, belongs to the ubiquitin-conjugating enzyme family. The family of ubiquitin-conjugating (E2) enzymes is characterized by the presence of a highly conserved ubiquitin-conjugating (UBC) domain. These domains accommodate the ATP-activated ubiquitin (Ub) or ubiquitin-like (UBL) protein via a covalently linked thioester onto its active-site residue. E2 enzymes act via selective protein-protein interactions with the E1 and E3 enzymes and connect activation to covalent modification. By doing so, E2s differentiate effects on downstream substrates, either with a single Ub/UBL molecule or as a chain. UBCH8 is highly similar in primary structure to the enzyme encoded by the UBE2L3 gene. It catalyzes the covalent attachment of ubiquitin or ISG15 to other proteins. UBCH8 functions in the E6/E6-AP-induced ubiquitination of p53/TP53 and promotes ubiquitination and subsequent proteasomal degradation of FLT3. At protein level, it is present in natural killer cells.

**配方:**PBS

**分子量:**19.6 kDa

**序列:**Met 1-Ser 153

**纯度:**> 95% by HPLC

**浓度:**

**内毒素:**<1.0 EU per 1 ug of protein (determined by LAL method)

**存储:** +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.