
本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human IGSF3 / EWI-3 Protein (His tag)

产品货号: TP06608

第三版

别名:EWI-3;LCDD;V8

描述:Immunoglobulin superfamily member 3 (IGSF3), also known as Glu-Trp-Ile EWI motif-containing protein 3, IgSF3 and EWI-3 and IGSF3, is a single-pass type I membrane protein which belongs to the EWI subfamily of the Ig Superfamily of molecules. IGSF3 has an overall structural similarity and strong sequence similarity to V7, a human leukocyte surface protein. IGSF3 contains eight Ig-like C2-type (immunoglobulin-like) domains. IGSF3 is expressed in a wide range of tissues with high expression in Placenta, kidney and lung. EWI-2 is part of a novel Ig subfamily that includes EWI-F (F2alpha receptor regulatory protein (FPRP), CD9P-1), EWI-3 (IGSF3), and EWI-101 (CD101). All four members of this Ig subfamily contain a Glu-Trp-Ile (EWI) motif not seen in other Ig proteins. Human IGSF3 shares 92 % aa identity with mouse IGSF3.

配方:PBS

分子量:126.6 kDa

序列:Met 1-Ala 1124

纯度:> 95% by HPLC

浓度:

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.