

---

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

---

Recombinant Human CNDP1 Protein (His Tag)

产品货号: TP06621

第三版

别名:CN1;CPGL2;HsT2308

**描述:**CNDP1, also known as carnosine dipeptidase 1, glutamate carboxypeptidase-like protein 2 (CPGL-2) or carnosinase 1 (CN1), is a member of the M20 metalloprotease family. The CNDP1 gene contains trinucleotide (CTG) repeat length polymorphism in the coding region, which has been demonstrated to be associated with susceptibility to developing diabetic nephropathy, for carnosine protection against the adverse effects of high glucose levels on renal cells. In humans, CNDP1 is secreted from the liver into the serum. In other mammals, including rodents, CNDP1 is expressed exclusively within the kidney and lacks a signal peptide. CNDP1 protein is a secreted homodimeric dipeptidase that specifically hydrolyzes L-carnosine (  $\gamma$ -alanyl-L-histidine), and is identified as human carnosinase expressed in the brain. CNDP1 has been associated with diabetic nephropathy in Europeans and European Americans, but not African-Americans. It was identified and confirmed as a risk factor, were cross-sectional and mostly in patients with type 2 diabetes. The polymorphisms of CNDP1 can be excluded as a risk factor for nephropathy in type 1 diabetes. In addition, CNDP1 is also suggested to be implicated in the actions of neuroprotection and neurotransmiting.

**配方:**PBS

**分子量:**55.3 kDa

**序列:**Ser27-His507

**纯度:**> 95% by HPLC

**浓度:**

**内毒素:**<1.0 EU per 1 ug of protein (determined by LAL method)

**存储:** +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.