

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human SPG21 Protein (GST tag)

产品货号: TP06623

第三版

别名: ACP33; BM-019; GL010; MAST

描述: Spastic paraplegia 21 (SPG21), also known as acid Cluster Protein 33 (ACP33) and Mast syndrome protein, is a member of the AB hydrolase superfamily. Human SPG21 is a 308 amino acid residue protein widely expressed in all tissues, including heart, brain, placenta, lung, liver, skeletal muscle, kidney and pancreas. SPG21 binds to the hydrophobic C-terminal amino acids of CD4 which are involved in repression of T cell activation via the noncatalytic alpha/beta hydrolase fold domain. SPG21 thus is proposed to play a role as a negative regulatory factor in CD4-dependent T-cell activation of CD4. Defects in SPG21 are the cause of spastic paraplegia autosomal recessive type 21, also known as Mast syndrome, a neurodegenerative disorder characterized by a slow, gradual, progressive weakness and spasticity of the lower limbs. Rate of progression and the severity of symptoms are quite variable. SPG21 is also associated with dementia and other central nervous system abnormalities.

配方: PBS

分子量: 61 kDa

序列: Met 1-Gln 308

纯度: > 95% by HPLC

浓度:

内毒素: <1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.