



本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human PF4V1 Protein (Fc Tag)

产品货号: TP06650

第三版

别名: CXCL4L1; CXCL4V1; PF4-ALT; PF4A; SCYB4V1

描述: The protein encoded by this gene is a chemokine that is highly similar to platelet factor 4. The encoded protein displays a strong antiangiogenic function and is regulated by chemokine (C-X-C motif) receptor 3. This protein also impairs tumor growth and can protect against blood-retinal barrier breakdown in diabetes patients.

配方: PBS

分子量: 34.6 kDa

序列: Met 1-Ser104

纯度: > 95% by HPLC

浓度:

内毒素: <1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.