



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human 14-3-3 eta / YWHAH Protein (GST tag)

产品货号: TP06669

第三版

别名:YWHA1

描述:YWHAH, the gene encoding the 14-3-3eta isoform, is highly expressed in retinal ganglion cells (RGC). YWHAH is a positional and functional candidate gene for both schizophrenia and bipolar disorder (BP). It is located on chromosome 22q12.3, a region which has been implicated by linkage studies in both BP and schizophrenia.

配方:PBS

分子量:55 kDa

序列:Gly 2-Asn 246

纯度:> 95% by HPLC

浓度:

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.