
本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human C6 / complement component 6 Protein (His tag)

产品货号: TP06689

第三版

别名:C6

描述:Neuromodulin, also known as Axonal membrane protein GAP-43, Growth-associated protein 43, Neural phosphoprotein B-50, pp46 and GAP43, is a cell membrane protein which belongs to the neuromodulin family. Neuromodulin / GAP43 contains one IQ domain. Neuromodulin / GAP43 is associated with nerve growth. It is a major component of the motile "growth cones" that form the tips of elongating axons. Neuromodulin / GAP43 is involved in neurite outgrowth, a crucial process for the differentiation of neurons. The sudden infant death syndrome (SIDS) is the main cause of postneonatal infant death and its cause is still unknown. Neuromodulin / GAP43 is a marker of synaptic plasticity and is critical for normal development of the serotonergic innervation. Neuromodulin / GAP43 is a major cortical cytoskeleton-associated and calmodulin binding protein that is widely and abundantly expressed during development, maintained in selected brain structures in the adult, and reinduced during nerve regeneration. CAP23 and GAP43 are functionally related intrinsic determinants of anatomical plasticity. These proteins function by locally promoting subplasmalemmal actin cytoskeleton accumulation.

配方:PBS

分子量:104 kDa

序列:Met 1-Ala 934

纯度:> 95% by HPLC

浓度:

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.