



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human CDK16 / PCTAIRE1 / PCTK1 Protein (GST tag)

产品货号: TP06703

第三版

别名:PCTAIRE;PCTAIRE1;PCTGAIRE;PCTK1

描述:The protein encoded by this gene belongs to the cdc2/cdkx subfamily of the ser/thr family of protein kinases. It may play a role in signal transduction cascades in terminally differentiated cells; in exocytosis; and in transport of secretory cargo from the endoplasmic reticulum. This gene is thought to escape X inactivation. Alternative splicing results in multiple transcript variants encoding different isoforms.

配方:PBS

分子量:82 kDa

序列:Met 1-Phe 496

纯度:> 95% by HPLC

浓度:

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.