

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Mouse KYNU / Kynureninase Protein (His tag)

产品货号: TP06708

第三版

别名:4432411A05Rik

描述:Genetic studies in mouse and human suggest that kynureninase activity may influence blood pressure and renal function. The gene coding kynureninase (KYNU) is also located on chromosome band 2q14-q23, where a linkage peak for essential hypertension was previously detected in the Chinese Han population. The results show that the rare KYNU variant Arg188Gln affects kynureninase activity and are consistent with the hypothesis that this mutation can predispose to essential hypertension.

配方:PBS

分子量:53.7 kDa

序列:Met 1-Ser 464

纯度:> 95% by HPLC

浓度:

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C 或者-70 ° C. 避免反复冻融.