



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human GKN1 / Gastrokine 1 Protein (His tag)

产品货号: TP06730

第三版

别名:AMP18;BRICD1;CA11;FOV;foveolin;UNQ489/PRO1005

描述:The protein encoded by this gene is found to be down-regulated in human gastric cancer tissue as compared to normal gastric mucosa.

配方:PBS

分子量:19.6 kDa

序列:Asn35-Asn199

纯度:> 95% by HPLC

浓度:

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.