
本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human APOA4 Protein (His Tag)

产品货号: TP06737

第三版

别名:

描述: Apolipoprotein is genetically associated with the risk of Alzheimer's disease (AD). The APOA1, APOC3, and APOA4 genes are closely linked and located on human chromosome 11. There was a decreased trend for levels of APOA1, APOC3, and APOA4 in AD patients. CONCLUSION: Low levels of APOA1, APOC3, and APOA4 are associated with risk of AD. APOA1, APOC3, and APOA4 should be developed as combined drugs for the therapy of AD. SNP(single nucleotide polymorphisms)in APOA1and APOA4 genes influences atherogenic characteristics of LDL particles in response to diet.

配方: PBS

分子量: 44.8 kDa

序列: Met 1-Ser396

纯度: > 95% by HPLC

浓度:

内毒素: <1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.