
本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Rat PNLIPRP1 / PLRP1 Protein (His Tag)

产品货号: TP06762

第三版

别名:Plrp1

描述:PNLIPRP1, also known as PLRP1, belongs to the AB hydrolase superfamily, Lipase family. PNLIPRP1 is structurally similar to PLRP2. However, these two proteins display different functional properties. PNLIPRP1 may function as inhibitor of dietary triglyceride digestion. It lacks detectable lipase activity towards triglycerides, diglycerides, phosphatidylcholine, galactolipids or cholesterol esters. PLRP2 hydrolyses milk fat with a lower catalytic efficiency than that of PL. PLRP2 activity, higher on homogenized than on native milk fat, is differently influenced by fatty acids and colipase depending on a proteolytic cleavage in the lid domain.

配方:PBS

分子量:52 kDa

序列:Met1-Val473

纯度:> 95% by HPLC

浓度:

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.