

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human IMP1 / IMPA1 Protein (His tag)

产品货号: TP06772

第三版

别名:IMP;IMPA

描述:The IMPA1 gene product is responsible for the final step of biotransformation of inositol triphosphate and diacylglycerol, two second messengers. Despite its many physiological functions, no clinical phenotype has been assigned to this gene dysfunction to date. Additionally, IMPA1 is the main target of lithium, a drug that is at the forefront of treatment for bipolar disorder.

配方:PBS

分子量:32.4 kDa

序列:Met 1-Asp277

纯度:> 95% by HPLC

浓度:

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.