

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human CAMK4 / CaMKIV Protein (GST tag)

产品货号: TP06781

第三版

别名:caMK;CaMK-GR;CaMKIV;IV

**描述:**Ca<sup>2+</sup>/calmodulin-dependent protein kinase 4 (CAMK ) belongs to the serine/threonine protein kinase family, and to the Ca<sup>2+</sup>/calmodulin-dependent protein kinase subfamily which is widely recognized as an essential enzyme implicated in the phosphoinositide amplification cascade. Ca<sup>2+</sup>/calmodulin dependent protein kinase (CAMK) can be activated by the intracellular increased Ca<sup>2+</sup> and then apt to combine with the target protein. Ca<sup>2+</sup>/calmodulin-dependent protein kinase 4 (CAMK ) is a multifunctional CaM-dependent kinase protein with limited tissue distribution, that has been implicated in transcriptional regulation in lymphocytes, neurons and male germ cells. All of the isoforms of this family, including myosin light chain kinase, phosphorylase kinase, CaMK1, CaMK and CaMK have EF-hand structure.

**配方:**PBS

**分子量:**79 kDa

**序列:**Met 1-Tyr 473

**纯度:**> 95% by HPLC

**浓度:**

**内毒素:**<1.0 EU per 1 ug of protein (determined by LAL method)

**存储:** +4 ° C 保存 (1-2 周). 长期保存在-20 ° C 或者-70 ° C. 避免反复冻融.