
本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human FAM20B / Gxk1 Protein (Fc tag)

产品货号: TP06798

第三版

别名:gxk1

描述:PSG6 is a pregnancy-specific glycoprotein (PSG). PSGs are secreted proteins which are produced by the rodent and primate placenta and play a critical role in pregnancy success. The levels of PSGs are highest during the third trimester of pregnancy, a time marked by the most profound suppression of MS disease attacks. PSGs regulate T-cell function. The regulation of T-cell function during pregnancy is likely the result of significant hormonal changes and may well involve immunoregulatory proteins derived from the placenta. Pregnancy specific glycoproteins (PSGs) are the most abundant placentally derived glycoproteins in the maternal serum. PSG1, PSG6, PSG6N, and PSG11 induce dose-dependent secretion of anti-inflammatory cytokines by human monocytes. Human and murine PSGs exhibit cross-species activity.

配方:PBS

分子量:71.4 kDa

序列:Ser31-Leu409

纯度:> 95% by HPLC

浓度:

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.