

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human Cbl-c / CBL-3 Protein (His & GST tag)

产品货号: TP06803

第三版

别名:CBL-3;CBL-SL;RNF57

描述:CBL proteins, such as Cbl-c, are phosphorylated upon activation of a variety of receptors that signal via protein tyrosine kinases. Through interactions with proteins containing SRC homology-2 (SH2) and SH3 domains, CBL proteins modulate downstream cell signaling. Cbl-c is a member of the Cbl family of E3 ubiquitin ligases. Expression of Cbl-c gene may be restricted to epithelial cells, and alternatively spliced transcript variants encoding multiple isoforms have been observed for Cbl-c gene.

配方:PBS

分子量:80.2 kDa

序列:Met 1-Ala474

纯度:> 95% by HPLC

浓度:

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.