

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human Galactolipase / PNLIPRP2 Protein (His tag)

产品货号: TP06806

第三版

别名:PLRP2

描述:Galactolipase, also known as PNLIPRP2, is a lipase with broad substrate specificity. It can hydrolyze both phospholipids and galactolipids. Galactolipase acts preferentially on monoglycerides, phospholipids and galactolipids. It also hydrolyses milk fat with a lower catalytic efficiency. The expressed galactolipase shows a lipolytic activity that is, however, only marginally dependent on the presence of colipase. The lipolytic activity of pancreatic extracts and human pancreatic juice on Labrasol is mainly due to the combined action of carboxyl ester hydrolase and galactolipase.

配方:PBS

分子量:51.5 kDa

序列:Met 1-Cys469

纯度:> 95% by HPLC

浓度:

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.