



本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human PTH2 / PTRH2 Protein (His tag)

产品货号: TP06823

第三版

别名:BIT1;CFAP37;CGI-147;IMNEPD;PTH;PTH 2;PTH2;PTRH2

描述:PTH2, also known as PTRH2, is a mitochondrial protein which belongs to the PTH2 family. PTH 2 is released during apoptosis from the mitochondria to the cytoplasm. The natural substrate for PTH 2 may be peptidyl-tRNAs which drop off the ribosome during protein synthesis. When in the cytoplasm, PTRH2 regulates the function of 2 transcriptional regulators, TLE5 and TLE1, thus promoting caspase-independent cell death.

配方:PBS

分子量:14.5 kDa

序列:Gly63-Tyr179

纯度:> 95% by HPLC

浓度:

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.