

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human LAP3 / Cytosol aminopeptidase Protein (His Tag)

产品货号: TP06828

第三版

别名:HEL-S-106;LAP;LAPEP;PEPS

描述:LAP3 (Leucine Aminopeptidase 3) is a Protein Coding gene. Diseases associated with LAP3 include Bacterial Vaginosis and Trichomoniasis. Among its related pathways are Arginine and proline metabolism and Metabolism. Gene Ontology (GO) annotations related to this gene include peptidase activity and aminopeptidase activity. An important paralog of this gene is NPEPL1.

配方:PBS

分子量:57.6 kDa

序列:Met 1-Ala519

纯度:> 95% by HPLC

浓度:

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C 或者-70 ° C. 避免反复冻融.