

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human ART3 Protein (His tag)

产品货号: TP06832

第三版

别名:ARTC3

描述:ART3 is an arginine-specific ADP-ribosyltransferase which belongs to the Arg-specific ADP-ribosyltransferase family. ART3 catalyzes a reversible reaction which modifies proteins by the addition or removal of ADP-ribose to an arginine residue to regulate the function of the modified protein. It is expressed specifically in testis. ART3 pseudogene is located on chromosome 11. ART3 was identified as a susceptibility gene for non-obstructive azoospermia (NOA). It is a novel therapeutic target in the treatment of NOA.

配方:PBS

分子量:39.6 kDa

序列:Met 1-Ser362

纯度:> 95% by HPLC

浓度:

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.