

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human Annexin A8 / ANXA8 Protein (His tag)

产品货号: TP06838

第三版

别名:ANX8

描述:DPEP2 (MBD-2) belongs to the membrane-bound dipeptidase family. There are three members of this family as membrane-bound dipeptidase-1 (MBD-1),membrane-bound dipeptidase-2 (MBD-2) and membrane-bound dipeptidase-3 (MBD-3).MBD-2 is expressed at highest levels in lung, heart, and testis and at some what lower levels in spleen.MBD-2 is membrane-bound through a glycosylphosphatidyl-inositol linkage and probably is a metalloprotease which hydrolyzes leukotriene D4 (LTD4) into leukotriene E4 (LTE4).It is generally recognized that rapid cleavage of LTD4 is important in inactivating the broncho-and vaso-constrictive effects of cysteinyl LTs in asthmatic and inflammatory processes.

配方:PBS

分子量:39.1 kDa

序列:Met 1-Pro327

纯度:> 95% by HPLC

浓度:

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.