



本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human CART / CARTPT Protein (Fc tag)

产品货号: TP06883

第三版

别名: CART

描述: This gene encodes a preproprotein that is proteolytically processed to generate multiple biologically active peptides. These peptides play a role in appetite, energy balance, maintenance of body weight, reward and addiction, and the stress response. Expression of a similar gene transcript in rodents is upregulated following administration of cocaine and amphetamine. Mutations in this gene are associated with susceptibility to obesity in humans.

配方: PBS

分子量: 36.3 kDa

序列: Met 1-Leu116

纯度: > 95% by HPLC

浓度:

内毒素: <1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.