

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human BOLA1 Protein (His tag)

产品货号: TP06886

第三版

别名: CGI-143

**描述:** Galactolipase, also known as PNLIPRP2, is a lipase with broad substrate specificity. It can hydrolyze both phospholipids and galactolipids. Galactolipase acts preferentially on monoglycerides, phospholipids and galactolipids. It also hydrolyses milk fat with a lower catalytic efficiency. The expressed galactolipase shows a lipolytic activity that is, however, only marginally dependent on the presence of colipase. The lipolytic activity of pancreatic extracts and human pancreatic juice on Labrasol is mainly due to the combined action of carboxyl ester hydrolase and galactolipase.

**配方:** PBS

**分子量:** 16.4 kDa

**序列:** Met 1-Pro137

**纯度:** > 95% by HPLC

**浓度:**

**内毒素:** <1.0 EU per 1 ug of protein (determined by LAL method)

**存储:** +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.