
本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human APLP1 / Amyloid-like protein 1 Protein (His tag)

产品货号: TP06893

第三版

别名: APLP

描述: APLP1, also known as amyloid-like protein 1, is a member of the highly conserved amyloid precursor protein gene family. APLP1 is a membrane-associated glycoprotein that is cleaved by secretases in a manner similar to amyloid beta A4 precursor protein cleavage. APLP1, together with APLP2, are important modulators of glucose. APLP1 may also play a role in synaptic maturation during cortical development. Alternatively spliced transcript variants encoding different isoforms have been described. APLP1 also is a mammalian homologue of amyloid precursor protein (APP). APP is a type I membrane protein that is genetically linked to Alzheimer's disease.

配方: PBS

分子量: 61.9 kDa

序列: Met 1-Glu580

纯度: > 95% by HPLC

浓度:

内毒素: <1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.