

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human FANCA / FACA Protein (His tag)

产品货号: TP06903

第三版

别名:FA;FA-H;FA1;FAA;FACA;FAH;FANCH

描述:FANCA is one of the six known Fanconi anemia gene products (FANCA, FANCC, FANCD2, FANCE, FANCF, and FANCG proteins). Fanconi anemia (FA) is a genetic disorder predisposing to aplastic anemia and cancer characterized by hypersensitivity to DNA-damaging agents and oxidative stress. FANCA associates with the I κ B kinase (IKK) signalsome via interaction with IKK2. Components of the FANCA complex undergo rapid, stimulus-dependent changes in phosphorylation, which are blocked by kinase-inactive IKK2.

配方:PBS

分子量:35.2 kDa

序列:Met 1-Cys 297

纯度:> 95% by HPLC

浓度:

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.