

---

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

---

Recombinant Human TRP1 / TYRP1 Protein (His tag)

产品货号: TP06910

---

第三版

别名:b-PROTEIN;CAS2;CATB;GP75;OCA3;TRP;TRP1;TYRP

**描述:** Tyrosinase-related protein 1, also known as TYRP1 or TRP1, is a melanosomal enzyme that belongs to the tyrosinase family and plays an important role in the melanin biosynthetic pathway. Mutations in this enzyme are the cause of rufous oculocutaneous albinism and oculocutaneous albinism type III. TYRP1 / TRP1 is involved in the oxidation of 5,6-dihydroxyindole-2-carboxylic acid (DHICA) into indole-5,6-quinone-2-carboxylic acid. This enzyme may regulate or influence the type of melanin synthesized. The expression of Tyrosinase-related protein 1 (TYRP1) is regulated by the microphthalmia-associated transcription factor (MITF). There is mounting evidence demonstrating that in addition to its role in eumelanin synthesis, TYRP1 is involved in maintaining stability of tyrosinase proliferation and melanocyte cell death.

**配方:** PBS

**分子量:** 52 kDa

**序列:** Met 1- Arg 471

**纯度:** > 95% by HPLC

**浓度:**

**内毒素:** <1.0 EU per 1 ug of protein (determined by LAL method)

**存储:** +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.